

All images used in this course are included in the course materials (see https://github.com/umcu-isi/programming-for-medical-imaging/releases/download/course-materials/course_materials.zip).

MicroDicom practicum (Windows only)

MicroDicom is a free tool to view medical images. These steps cover installation, basic tools, and measurements.

1. Make sure you have installed MicroDicom on your computer.
2. 2D image manipulation. Click File and select Open (or drag a file) to load the DICOM file "Knee.dcm", which is indicated as Anonymous.
 - Look for the most regular functions to manipulate images:
 - Click the Zoom icon or hold the right mouse button to make the image bigger or smaller.
 - Select the Pan icon (hand symbol) to move the image within the viewport.
 - Adjust the brightness and contrast by clicking the Window/Level tool and dragging your mouse over the picture.
 - Create some annotations (lines, curves, circles, etc) and delete them again.
 - Find out how to measure the angle of 2 lines.
 - Apply different color scales, display the knee for instance in the heated object scale (HOT).
 - The image contains a metal wedge (aluminium). What could be a reason to add this device in an X-ray?
 - Each step of the wedge is 1cm high, 4cm width. Find the measurement tool and measure the length of the wedge, e.g. of 10 steps (10cm) height. With this knowledge, you should be able to calculate the pixel size of this X-ray. Show the calculation to a teaching assistant (TA), or ask for help.
 - Find the stored pixel size in the DICOM header. Compare this with the results of your distance measurement(s). What could be causes for the differences?
3. Next, drag a full stack of CT slices to MicroDicom. Choose any set you want.
 - Use your mouse wheel to scroll up and down through the image slices. Check the effect and differences if you use the arrow buttons and the PgUp PgDn buttons to scroll through the images.
 - Use the Angle tool to draw lines and find the angle between them.
 - Place the ROI (Region Of Interest) tool on a specific spot to see average pixel values.
 - Choose View->Dicom Tags Pane to show the DICOM header information of the current image on the right.
 - Scroll through the 3D image and try to find tags that change if you scroll.
 - You can add some tags to Favorites with your mouse (right mouse button click).
 - Open the CT pelvis and add the tags "(0020,1041)", "(0020,0032)", and "(0018,0050)" to Favorites. Make the Favorites Tags active (right below) and start scrolling again. Try to explain the changes in the values and ask for help if you don't understand it.

Osirix practicum (macOS only)

Osirix Lite is a free tool to view and study medical images. These steps cover installation, basic tools, and measurements.

1. Make sure you have installed Osirix Lite on your computer.
2. 2D image manipulation. Use Import to add images to your database. Open the DICOM file "Knee.dcm", which is indicated as Anonymous.
 - Look for the most regular functions to manipulate images:
 - Hold the right mouse button to make the image bigger or smaller.
 - Click the Command button and the left mouse button to move the image within the viewport (panning).
 - Adjust the brightness and contrast by using the middle mouse button.
 - Click on "2D Viewer"; this will switch to a viewing mode with options.
 - Play with the mouse button functions of the top row.
 - Create some annotations (lines, curves, polygons, etc) and delete them again.
 - Find out how to measure the angle of 2 lines.
 - Apply different color scales, display the knee for instance in a heated object scale (HotMetal).

- The image contains a metal wedge (aluminium). What could be a reason to add this device in an X-ray?
 - Each step of the wedge is 1cm high, 4cm width. Find the measurement tool and measure the length of the wedge, e.g. of 10 steps (10cm) height. With this knowledge, you should be able to calculate the pixel size of this X-ray. Show the calculation to a teaching assistant (TA), or ask for help.
 - Look for the Meta-Data of this image and find the stored pixel size in the DICOM header. Compare this with the results of your distance measurement(s). What could be causes for the differences?
3. Next, open the full stack of CT slices "James Bond" in the database.
- Use your mouse wheel to scroll up and down through the image slices.
 - Go to the 2D Viewer mode again.
 - Make all image information visible in the viewer.
 - Change the viewing direction (sagittal, coronal, axial) and scroll through the slices.
 - Check some pixel values (air, bone, etc) by moving your mouse over the image.
 - Place the ROI (Region Of Interest) tool on a specific spot to see some statistical values.
 - Choose File->Meta-Data to open the Tags Pane. This way, the DICOM header information of the current image will be shown on the right.
 - Scroll through the 3D image and try to find tags that change if you scroll. Note: you may have to enlarge the DICOM pane to see the actual values.
 - Use the Search field to look for "ImagePositionPatient" and try to understand what you see when you scroll through the image. Ask for help if you don't understand it.

ITK-SNAP practicum

ITK-SNAP is also a free tool to view and study medical images. This application can also deal with images in Nifti format, which we will use throughout the PMI course as well.

1. Make sure you have installed ITK-SNAP on your computer.
2. Click File and select Open (or drag a file) to load the Nifti file "ct_jaw_slice.nii.gz".
3. Look for the most regular functions to manipulate images:
 - Use the right mouse button to zoom in and out.
 - Panning can be done with the left mouse button.
 - Try to find the tool that handles the contrast (window/level). Adjust the brightness and contrast and play with the available options.
4. Check the Metadata of the image. Do you see the differences with the header information of DICOM images?
5. Try to find the Annotation tool, draw some lines and look at the information shown in the image.
 - Change the color of the lines.
 - Delete an annotation.
6. Load image "ct_jaw.nii.gz" as new Main image.
7. Scroll in all 3 visible windows (viewports) and make sure you understand what is being displayed in every viewport as you scroll (use mouse wheel and arrow buttons).
8. Load "segment_1.nii.gz" as Segmentation (drag this file into ITK-SNAP, a window will pop-up).
 - Scroll through all slices and try to find the abnormality of this patient.
9. Unload the mask as segmentation, and create your own segmentation(s) using the Paint brush mode.
 - Use multiple colors for different labels.
 - Try to calculate the volume of your own drawings.

Online web viewer

You can also use a simple web viewer to view Nifti and/or DICOM images.

- Go to <https://mix.isi.uu.nl/viewer/> and drag the image of your choice into the viewer.
- Play around with your mouse to manipulate the view.